ASSUT EUROPE

SUTURE CHIRURGICHE

More than surgery

**PRODOTTI SPECIALI
SPECIAL PRODUCTS
PRODUCTOS ESPECIALES**

LAPAROMESH® II

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2- Cellulosa.

 **LAPAROMESH® II** is a new semi-absorbable composite design mesh made from a network of a synthetic absorbable monofilament in poly (glycolide-co- Σ -caprolactone), polypropylene monofilament and a cellulose film on the visceral side. Laparomesh is designed for the prevention of incisional hernias in patients with compromised wound healing. Healing is improved as the growth of connective tissue is promoted, increasing the width and height of the linea alba stronger tensile strength over that recorded for closures achieved through simple suture. Laparomesh should be used only in closure procedures performed in clean and/or clean - contaminated surgery.

STRUCTURE

1- Vertical arm made from synthetic absorbable monofilament in poly (glycolide-co- Σ -caprolactone) and polypropylene.

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


SPECIFICATIONS

Initial density of the mesh: 85-90 g/m²
 Pore diameter: 2,0 ± 0,3 mm
 Thickness: 500 ± 50 µm
 Content of Polypropylene: 30 g/m²

ABSORPTION

Complete absorption of poly (glicolide-co-ε-caprolattone) between 90 and 110 days after implantation only a lightweight mesh remains

INDICATIONS

Laparomesh® is used in patients with high risk factors:

Oncological patients
 Diabetics
 Immunodepressed patients
 Patients with morbid or extreme obesity
 Elderly patients
 And in the following specialties:
 Vascular Surgery
 General Surgery
 Urological Surgery
 Gynaecological Surgery

BENEFITS

Excellent flexibility
 Excellent handling
 Reduces the possibility of inflammation
 The large-porous mesh enhances the formation of connective tissue

METHOD OF USE

Laparomesh® II is implanted at the midline by including its vertical component between the edges of the aponeurosis. Fascial closure (suture+mesh) is achieved with a running polydioxanone or polypropylene suture (monofilament Assuofil, Assupro) or Filbloc.
 The horizontal mesh component may be placed, or not, in contact with the visceral peritoneum since adhesions are avoided by its cellulose coating.

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